mini-GUIDE



USING SPARKvue SOFTWARE WITH PASCO WEATHER SENSOR TO GATHER ONLY TEMPERATURE DATA (OR ANY OTHER SINGLE MEASUREMENT OF CHOICE)

STEP 1:

Connect to Pasco Weather Sensor:

See Mini-Guide called "Set Up the Pasco Weather Sensor" to show you how to connect SPARKvue software to the Pasco Weather Sensor by Bluetooth. Click 'Build New Experiment' as instructed.

STEP 2:

Select the single window display type.



STEP 3:

Select the line graph to display your data.

STEP 4:

Click 'Select Measure' for Y-axis and choose 'Temperature'.

The X-Axis measure defaults to 'Time', which is what we want for this experiment.

Choose how often you want to record a measure (2 hz = 2 times per second)



Information 🗟 🕤 At A Glance





STEP 5:

Click the 'tools' button and 'data properties' to select whether temperature is measured according to Celsius (C), Farenheit (F) or Kelvin (K) scales.

]	1: Untitled		8 6
-	Experiment Too	ls	?	
Temperature (10)		Menage Runs		
0	- [Manage Images Galculated Data		
-6		Data Properties		
-10	-		ż	8 9
<u>k.</u> *		Done		Palative Humidity (%)
Periodic: 2 Hz	% • Temperature 22.9	Start 00:00:00.0		

STEP 6:

Click 'Start' to begin recording data. The same button becomes the Stop button to let you stop data collecting.¹



¹ All content comes from <u>https://www.pasco.com/products/sensors/environmental/ps-3209</u> <u>https://www.pasco.com/products/software/sparkvue</u>

Information 🗟 🗿 At A Glance